

CODE: WC740
ISSUED: 2009 NOV
ABLE RATING **AMBER**



★ ALTERNATIVES OR SIMILAR ★

INFORMATION SHEET

LIQUID LAB CLEANER

DESCRIPTION:

LIQUID LAB CLEANER is a specially designed to clean laboratory glassware and slides. Formulated for easy rinsing, leaving glassware streak free.

RECOMMENDED USE:

1 1/2ml per litre of water at 45-55°C. Thoroughly wash glassware, using a brush if necessary. Rinse in clean water and allow to air dry.


SAFETY PRECAUTIONS:

WARNING: CORROSIVE - Avoid contact with skin and eyes.

FIRST AID:

PACKAGING:

Liquid Lab Cleaner is available in the following pack sizes



CONDITIONS OF SALE & USAGE (STANDARD DISCLAIMER)

This product should not be used for any purpose or in any manner contrary to the labelling, SDS, MSDS or Information Sheet unless authorised to do so under appropriate legislation. By purchasing and using this product the purchaser assumes responsibility for any and all results that may occur and shall not hold the manufacturer liable for any loss or injury however occurring.

BENEFITS & TIPS

- DO NOT MIX WITH OTHER CHEMICALS
- VERY HIGHLY CONCENTRATED.



PROUDLY AUSTRALIAN OWNED & MANUFACTURED BY ABLE WESTCHEM
ABN 82 009 353 182. 47 BANKSIA RD, WELSHPOOLWA 6106
Ph 08 9458 1604 Fax 08 9458 1222 Email: info@westchem.com.au
Web 1: www.westchem.com.au Web 2: www.envirocareearth.com.au

ALWAYS READ MSDS

